

# DIGITAL BLOOD PRESSURE MONITOR

## UA-774

**IHB**



### Dual Memory Monitor

- Irregular Heart Beat (IHB) Indicator
- Blood Pressure Classification Indicator
- Dual 30 Memory + Average Reading
- SlimFit Comfortable Cuff
- Illuminated Start Button
- Standard AC Adapter Inlet



**A&D**  
A&D Medical  
*Technology You Can Trust*  
<http://www.aandd.co.jp>

# DIGITAL BLOOD PRESSURE MONITOR UA-774

## Irregular Heart Beat (IHB) Indicator

This lights up if IHB is detected during the measurement. Helps users to find an irregular heartbeat, which may cause cardiovascular disease. (Note: This is for warning purposes only. A&D recommends contacting a physician if the user sees this indicator frequently)

## Blood Pressure Classification Indicator

World Health Organization (WHO) classifies blood pressure ranges into 6 stages. This monitor indicates your position (Optimal to Severe Hypertension) automatically after each measurement is complete.

## Dual 30 Memory + Average Reading

The UA-774 is equipped with dual 30-memory capability. Patients can track vital signs data for almost a month (one measurement per day). Up to 30 sets of data are automatically stored into the built-in memory of the unit. Averaged value of stored data will be shown for quick information and better appreciation of BP trend.

## SlimFit Comfortable Cuff

A new innovative edge free internal bladder "EFB" promises even pressure on the arm during the measurement. This provides the user with a higher level of comfort while the arm is under pressure.

## Illuminated Start Button


A&D has included an innovative idea for ease-of-use for two persons. The UA-774 has illuminated buttons and the button (1 or 2) is lit up while a measurement is in progress.

## Standard AC Adapter Inlet

The monitor is equipped with an AC inlet for external power source. You may be free from worrying about battery condition and enjoy economical use of the monitor for your daily use. (AC adapter is an option and not included in the package.)

## WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart below.



**A&D**  
A&D Medical  
*Technology You Can Trust*  
Tokyo JAPAN  
<http://www.aandd.co.jp>

\*UA774CA-ADCC-01-PR1-03a08  
REF\*CARD-UA774